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Bacteriology Lab 1

By: **Dr. A. Mohammadi**

Department of Biology,
Faculty of science,
University of Alzahra

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Isolation and Identification of Streptococci

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



Alpha Hemolysis



Gamma Hemolysis

■ **Beta hemolysis** refers to a clear, colorless zone surrounding the colony, where a complete lysis of the red blood cells by the hemolysins has occurred.


■ **Alpha hemolysis** appears as a zone of partial hemolysis surrounding the colony, often accompanied by a **greenish** or **brown** discoloration due to oxidation of hemoglobin to met-hemoglobin by hydrogen peroxide

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■ **Gamma hemolysis** refers to no hemolysis or discoloration of the agar surrounding the colony.

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
Bacitracin Sensitivity Test:


Results:

Positive test: any zone of inhibition around the disc.

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

4-5 BILE ESCULIN TEST RESULTS ♦ This plate was inoculated with two Gram-positive cocci. The organism on the bottom is bile esculin-positive. The organism on the top is negative.

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TABLE OF RESULTS

| Result | Interpretation | Symbol |
|---|---|--------|
| Medium is darkened within 48 hours | Presumptive identification as a member of Group D <i>Streptococcus</i> or <i>Enterococcus</i> | + |
| No darkening of the medium after 48 hours | Presumptive determination as not a member of Group D <i>Streptococcus</i> or <i>Enterococcus</i> | - |

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Pyrrolidonyl-beta-naphthylamide (PYR) test

S.pyogenes possesses the enzyme **L-pyroglutamyl aminopeptidase** which hydrolyzes **Pyrrolidonyl-beta-naphthylamide** with the formation of free β - naphthylamine, which combines with cinnamaldehyde reagent to form a bright- red end product.


Method :

Rub a small amount of a colony to be tested on a surface of a filter paper impregnated with PYR (available commercially).

Add one drop of freshly reconstituted detector reagent N,N- dimethylaminocinnamaldehyde, with detergent (available commercially), and observe for a red color within 5 minutes.

Hippurate hydrolysis test:

Group B streptococci contain the enzyme hippuricase, which can hydrolyze hippuric acid. Other β streptococci lack this enzyme.

The production of hippuricase results in the hydrolysis of sodium hippurate with the formation of sodium benzoate and glycine.

The hydrolysis of hippurate can be detected by one of two methods:

- ❑ *Standard test for benzoate using ferric chloride*
- ❑ *Rapid test for glycine using ninhydrin reagent*

Standard test for benzoate using ferric chloride

Media and reagent:

Sodium hippurate medium

Brain Heart infusion broth.

1% Sodium hippurate.

Ferric chloride reagent

FeCl_3 12g.

2% aqueous HCl 100 ml.

Procedure :

Incubate a tube of sodium hippurate medium with the organism.


Incubate for 20 hours or longer at 37 °C.

Centrifuge the medium to pack the cells and pipette 0.8 ml of the clear supernatant in to a clean test tube.

Add 0.2 ml of the ferric chloride reagent and mix.

Result:

A heavy precipitate will form upon addition of FeCl_3 reagent.

Rapid test for glycine using ninhydrin reagent

Media and reagent:

Sodium hippurate medium

1% Sodium hippurate.
Distilled water.

Ninhydrin reagent

Ninhydrin 3.5 g.
1: 1 Acetone/butanol 100 ml.

Procedure :

Incubate a large loop of a β streptococci culture in to 0.4 ml of the sodium hippurate medium and incubate for 2 hours at 37 °C.
Add 0.2 ml of the ninhydrin reagent.

Result:

The presence of a deep purple color after the addition of ninhydrin reagent within 10 minutes indicate the presence of glycine.

Quellung Reaction

Around the turn of the century (1900) German physician and bacteriologist Friedrich Neufeld discovered that antibodies against specific pneumococcal capsular antigens could be produced and used for typing *Streptococcus pneumoniae* organisms.

These antibodies, when bound to the cell wall antigen produced a clearing around each individual cell with the appearance of the capsule having swollen. The word "**Quellung**" is German for '**swelling**' however this is a misnomer as the capsule does not swell but simply appears enlarged with the clear zone produced by the bound antibodies.


The clearing is best visualized by using a stain to enhance contrast between the clear zone of bound antibody and the surrounding material.

Procedure:


- ⊕ Colony isolate, is placed on a clean microscope slide and allowed to dry.
- ⊕ A drop of *Polyvalent Quellung antisera is applied to the slide and mixed with a drop of methylene blue stain.
- ⊕ A coverslip is placed onto the mixture and allowed to incubate/react at room temperature for 15-20 minutes.
- ⊕ The slide is then examined under the microscope for distinct clearing around the cells.
- ⊕ One may have to search for a microscopic field which best shows the clearing. A distinct zone of clearing indicates a positive Quellung reaction and confirms the identification as *Streptococcus pneumoniae*.

***Polyvalent Antisera**- contains all known serotypes to identify any pneumococcus. Specific or individual antisera can be employed for scientific study of individual serotypes.

Quellung Reaction

Fluorescent antibody (FA) staining


***Streptococcus pneumoniae*, FA stain showing its antiphagocytic capsule (CDC).**

***Streptococcus faecalis* broth (SF broth) :**

Media and reagent:

Tryptone
D-glucose
Monopotassium phosphate
Dipotassium phosphate
Sodium azide
Sodium chloride
Bromocresol purple
pH= 6.9

Principle:

The tryptone and dextrose supply the nutrients required for the growth of enterococci. Sodium azide exhibits a bacteriostatic effect on gram-negative bacteria through its inhibitory action on enzymes in the electron transport system. Bromocresol purple serves as a pH indicator.

Pyrrolidonyl-beta-naphthylamide (PYR) test

Enterococcus spp. Posses the enzyme **L-pyroglutamyl aminopeptidase** which hydrolyzes **Pyrrolidonyl-beta-naphthylamide** with the formation of free β - naphthylamine, which combines with cinnamaldehyde reagent to form a bright- red end product.

Method :

Rub a small amount of a colony to be tested on a surface of a filter paper impregnated with PYR (available commercially).

Add one drop of freshly reconstituted detector reagent N,N- dimethylaminocinnamaldehyde, with detergent (available commercially), and observe for a red color within 5 minutes.

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α -hemolytic Streptococci

Optochin Sensitivity Test:

Definitive test to differentiate between

- *S. pneumoniae* & viridans Streptococci


Principle:

S.pneumoniae is inhibited by less than 5 $\mu\text{g/ml}$ Optochin reagent (ethylhydroxycupreine hydrochloride) giving a zone of inhibition.


Optochin Sensitivity Test:

Procedure:

1. Inoculate blood agar plate with the test organism.
2. Aseptically apply Optochin disc onto the center of the streaked area.
3. Incubate the plate at 37°C for 24 hrs.
4. Accurately measure the diameter of the inhibition zone around the disc.

Optochin Sensitivity Test:

Results:

Sensitive with inhibition zone equal or more than 14 mm when disk with diameter equal 6 mm or 16 when 10 mm disk is used.

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