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# Bacteriology Lab 1

By: **Dr. A. Mohammadi**

Department of Biology,  
Faculty of science,  
University of Alzahra

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# Isolation Of Enterobacteriaceae

- ❖ To isolate Enterobacteriaceae and Pseudomonas, specimens from the infected site are plated out on any one of a large number of selective and differential media such as EMB agar, Endo agar, Deoxycholate agar, MacConkey agar, Hektoen Enteric agar, and XLD agar.
- ❖ **Xylose Lysine Deoxycholate** agar is selective for gram-negative bacteria XLD agar contains: sodium desoxycholate ,sugars lactose and sucrose, the amino acid L-lysine, sodium thiosulfate, ferric ammonium citrate and the pH indicator phenol red.
- ❖ sodium desoxycholate, which inhibits the growth of gram positive bacteria but permits the growth of gram-negatives.

- ❖ If the gram-negative bacterium ferments lactose and/or sucrose, acid end products will be produced and cause the colonies and the phenol red in the agar around the colonies to turn yellow. If lactose and sucrose are not fermented but the amino acid lysine is decarboxylated, ammonia, an alkaline end product will cause the phenol red in the agar around the colonies to turn a deeper red. Sometimes the sugars are fermented producing acid end products and lysine is broken down producing alkaline end products in this case some of the colonies and part of the agar turns yellow and some of the colonies and part of the agar turns a deeper red.
- ❖ If hydrogen sulfide is produced from thiosulfate reduction, part or the entire colony will appear black.

✓ Typical reactions for some of the Enterobacteriaceae and Pseudomonas are:

- ❖ *Escherichia coli*: flat yellow colonies; some strains may be inhibited.
- ❖ *Enterobacter* and *Klebsiella*: mucoid yellow colonies.
- ❖ *Proteus*: red to yellow colonies; may have black centers.
- ❖ *Salmonella*: usually red colonies with black centers.
- ❖ *Shigella* and *Pseudomonas*: red colonies without black centers.

# Isolation Of *Pseudomonas*

- ❖ **Pseudosel agar** is selective for *Pseudomonas aeruginosa* and also stimulates *P. aeruginosa* to produce its characteristic pigment as well as fluorescent products. *Pseudomonas aeruginosa* will typically produce a green to blue water-soluble pigment on this agar and will also fluoresce when the plate is placed under a short wavelength ultraviolet light.
- ❖ Pseudosel agar contains cetrимide which inhibits most bacteria other than *Pseudomonas aeruginosa*. It also contains potassium sulfate and magnesium chloride which stimulate *P. aeruginosa* to produce the pyocyanin (a blue , non fluorescent pigment soluble in water and chloroform)
- ❖ Pseudosel agar can be used to stimulate the production of pigment and fluorescent products.

- ✓ **Pigment production on Pseudosel agar**
  - ❖ *Pseudomonas aeruginosa* produces a green to blue, water soluble pigment called pyocyanin.
  
- ✓ **Fluorescence under ultraviolet light on Pseudosel agar**
  - ❖ It produces a product called fluorescein that will fluoresce under short wavelength (254nm) ultraviolet light.
  
- ✓ **Odor**
  - ❖ Most of the Enterobacteriaceae have a rather foul smell; *Pseudomonas aeruginosa* produces a characteristic fruity or grape juice-like aroma due to production of an aromatic compound called aminoacetophenone.
  
- ✓ **Fermentation of glucose.**
  - ❖ All of the Enterobacteriaceae ferment the sugar glucose; *Pseudomonas aeruginosa* and other non fermentative gram negative rods will not.

**Pigment production on Pseudoseal agar**  
**Fluorescence under ultraviolet light on Pseudoseal agar**

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


# Oxidase test


16

**5-21** OXIDASE TEST ON BACTERIAL GROWTH ♦ A few drops of reagent on oxidase-positive bacteria will produce a purple-blue color immediately.

**5-22** OXIDASE SLIDE TEST ♦ Positive results with this test should appear within 20 seconds. The dark blue is a positive result. No color change is a negative result. The organism in the top right panel demonstrates a negative result (not blue) but appears yellow because it is pigmented. (BBL™ DrySlide™ systems available from Becton Dickinson, Sparks, MD.)


ASM MicrobeLibrary © Cathcart, Kramer, Shields

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


**5-25** INCUBATED NITRATE BROTH BEFORE THE ADDITION OF REAGENTS ♦ Numbered 1 through 5 from left to right, these are Nitrate Broths immediately after incubation before the addition of reagents. Tube 2 is an uninoculated control used for color comparison. Note the gas produced by tube 5. It is a known nonfermenter and, therefore, will receive no reagents. The gas produced is an indication of denitrification and a positive result. Tubes 1 through 4 will receive reagents. See Figure 5-26.

26


### 5-26 INCUBATED NITRATE BROTH AFTER ADDITION OF REAGENTS ♦

After the addition of reagents, tube 1 shows a positive result. Tube 3 and tube 4 are inconclusive because they show no color change. Zinc dust must be added to tubes 2 (control), 3, and 4 to verify the presence or absence of nitrate. See Figure 5-27.

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**5-27** INCUBATED NITRATE BROTH AFTER ADDITION OF REAGENTS AND ZINC ♦ Finally, a pinch of zinc is added to tubes 2, 3, and 4 because they have been colorless up to this point. Tube 2 (the control tube) and tube 3 turned red. This is a negative result because it indicates that nitrate is still present in the tube. Tube 4 did not change color, which indicates that the nitrate was reduced by the organism beyond nitrite to some other nitrogenous compound. This is a positive result.

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2. Incubate at  $35 \pm 2^\circ\text{C}$  for 24 hours

3. Examine the Durham tube

If there is a bubble and organism is NOT a fermenter, organism reduces nitrate to  $\text{N}_2$ . Test is completed.

4. If no bubble, add 8 drops Reagent A and 8 drops Reagent B—Mix well.

If red, organism reduces nitrate to nitrite. Test is completed.

5. If not red, add zinc.

If red, organism does not reduce nitrate. Test is completed.

If not red, organism reduces nitrate to something other than nitrite. Test is completed.

### TABLE OF RESULTS

Result	Interpretation	Symbol
Gas (Nonfermenter)	Denitrification—production of nitrogen gas ( $\text{NO}_3^- \rightarrow \text{NO}_2^- \rightarrow \text{N}_2$ )	+
Gas (Fermenter, or status is unknown)	Source of gas is unknown; requires addition of reagents	
Red color (after addition of reagents A and B)	Nitrate reduction to nitrite ( $\text{NO}_3^- \rightarrow \text{NO}_2^-$ )	+
No color (after the addition of reagents)	Incomplete test; requires the addition of zinc dust	
No color change (after addition of zinc)	Nitrate reduction to nongaseous nitrogenous compounds ( $\text{NO}_3^- \rightarrow \text{NO}_2^- \rightarrow$ nongaseous nitrogenous products)	+
Red color (after addition of zinc dust)	No nitrate reduction	-

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
TABLE 4-5 MacConkey Agar Results and Interpretations

## T A B L E O F R E S U L T S

Result	Interpretation	Presumptive ID
Poor growth or no growth (P)	Organism is inhibited by crystal violet and/or bile	Gram-positive
Good growth (G)	Organism is not inhibited by crystal violet or bile	Gram-negative
Pink to red growth with or without bile precipitate (R)	Organism produces acid from lactose fermentation	Probable coliform
Growth is "colorless" (not red or pink) (C)	Organism does not ferment lactose. No reaction	Noncoliform

**4-9** **MACCONKEY AGAR** + This MacConkey Agar was inoculated with four members of *Enterobacteriaceae*, all of which show abundant growth. The top organism and the one on the right are lactose fermenters, as evidenced by the pink color. Note the bile precipitate around the top organism. The organisms on the bottom and on the left produced no color, so they do not appear to be lactose fermenters.

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